

Jaco UltraLite Storage

Equip your carts for bedside medication administration. Jaco UltraLite Storage is lightweight, durable, and secure – and meets the most stringent medication delivery requirements. Mix and match components to create the storage solution you need.

Single-Drawer Cabinet

15.5"w x 14"d x 7"h



3" drawer
Made from the same lightweight aluminum as our cart bodies



3 Med Bins
Optional transparent covers secure bin contents, and provide visibility and easy access

Double-Drawer Cabinet

15.5"w x 14"d x 11"h



2 3" drawers



6 Med bins



Drawer & Med bins



7" drawer

Drawer Divider Kit

Divide drawers into 2-9 compartments



3 Compartments



6 Compartments



9 Compartments

Locking Systems

Push-Button Lock:

Simple, effective security for drawers

Touch-Pad Lock:

4-digit lock for drawers and med bins – includes on-board keypad

Keyboard E-Lock:

For drawers and med bins. Enter PIN via computer keyboard.
Track access history via Jaco FleetManager software



Need a custom storage solution? Call us.

Let our engineers help you create a storage solution to support your unique workflow. **Call 800.649.2278 to get started.**