

New Jaco One SL Storage System


Equip Jaco One for bedside medication administration – with our lightest, slimmest, most flexible storage system. New Jaco One SL Storage is lightweight, durable, secure – and engineered to meet the most stringent medication delivery requirements. Mix and match components to create the storage configuration you need.

Single-Drawer Cabinet

13.2" w x 13" d x 4.3" h



3" drawer



3 3.5" w Med Bins




7.2" w & 3.5" w Med Bins

Optional transparent Med Bin lids secure contents while providing easy access and visibility.

Double-Drawer Cabinet


13.2" w x 13" d x 8.5" h



2 3" drawers



Med Bins



Drawer & Med Bins



7" drawer

Triple-Drawer Cabinet

13.2" w x 13" d x 12.75" h



3 3" drawers



Med Bins



Drawers & Med Bins

Drawer Divider Kit

Divide Drawers into 2-9 Compartments



3 Compartments



6 Compartments



9 Compartments

Locking System

Push-Button Lock:

Simple, effective mechanical security for drawers

Touch-Pad Lock:

4-digit electronic lock for drawers or med bins – includes on-board keypad



Need a custom storage solution? We're ready to help.

Let our engineers create a storage solution to support your specific workflow requirements. Call 800.649.2278 to get started.