



NOW PART OF
GCX
Healthcare Solutions



EVO Documentation Carts

Fully-equipped mobile EMR workstations help nurses focus on accurate documentation and better care.

Lightweight and maneuverable

Up to 45% lighter and 29% easier to roll, with a compact profile that's easy to maneuver through tight spaces and crowded patient rooms.

Research-based ergonomics

Instant, silent, one-handed Comfort-Glide height adjustment; easy-to-access side-opening EVO ClearView Storage drawers; generous worksurface and well-positioned tools.

Superior infection control; UV-ready

The high-gloss, bright white antimicrobial finish always looks fresh, wipes clean easily and stands up to UV-B and UV-C sterilization.

Versatile and adaptable

Equip carts for any documentation, rounding or point-of-care workflow—and modify easily when requirements change.

A smart investment

EVO documentation carts built to last from aircraft-grade metal, and backed by our industry-leading 10-year warranty.

Popular EVO Documentation Cart Models

Choose below or build your own. All EVO documentation workstations are fully customizable to meet the demands of any workflow.

LCD Documentation Cart

Full-sized, full-shift powered medical cart for virtually any documentation workflow.

- EVO LCD cart with 500 Wh LiFePO₄ Power System*
- ClearView 7" drawer cabinet (e-Lock)
- Corded scanner holder
- Label/Wristband printer mount
- 5" diameter wipes holder

**Also available with POWER BLADE hot-swappable LiFePO₄ 476 Wh power system*



Locking Laptop Documentation Cart

Lightweight, durable powered laptop documentation cart

- EVO locking laptop cart with 250 Wh LiFePO₄ power system (EVO-10-L250)
- ClearView 3"/3" drawer cabinet (eLock)
- Wireless barcode scanner base mount
- Label/wristband printer mount
- 5" diameter wipes holder

Podium Laptop Documentation Cart

Sturdy, ultra-lightweight documentation cart —an easy-to-share, cost-effective solution

- EVO podium laptop cart (EVO-PODIUM)
- Laptop bar lock/security bridle
- ClearView 3" drawer cabinet (non-locking)
- Corded barcode scanner holder
- 5" diameter wipes holder



Add accessories and storage to create your own custom solution

EVO cart accessories

Wipes holders	Wire baskets	Wireless scanner base mount	Corded scanner holder	Label/wristband printer mounts
				
Worksurface utility cups	Medical gloves box holders	Height-adjustable VESA mount	Laptop bar lock/security bridle	Dual-wheel casters
				

EVO ClearView Storage options


Single Drawer

14.5"w x 9.25"d x 4.5"h




Double Drawer

14.5"w x 9.25"d x 8.25"h



Triple Drawer

14.5"w x 9.25"d x 11.9"h



Single Drawer

12.5"w x 7.5"d x 3.5"h




1 Large Bin

10"w x 7"d x 3.2"h



2 Long Bins

10.5"w x 3"d x 3"h




3 Mini Bins

3.5"w x 7.3"d x 2.75"h



Deep Drawer

12.5"w x 7.5"d x 7.1"h



* All drawer dimensions are interior dimensions.



NOW PART OF
GCTX
Healthcare Solutions

140 Constitution Blvd. | Franklin, MA 02038
800.649.2278 | sales@gcx.com