



NOW PART OF  
**GCTX**  
Healthcare Solutions



# Enable healthcare anywhere with Jaco Telehealth solutions

EVO telemedicine carts and wall mounted telemedicine workstations let clinicians deliver state-of-the-art telemedicine in any care setting.

## **Superior patient interaction**

Slim, compact, lightweight carts and flexible, fully adjustable wall mounts enable optimal real-time interaction between patients and providers.

## **Flexible and adaptable**

Designed to support best-of-breed videoconferencing and telemedicine equipment today, and new technologies as they emerge.

## **State-of-the-art infection control**

EVO's high-gloss, non-porous antimicrobial powder coat finish wipes clean easily and stands up to UV-B/UV-C sterilization.

## **Durable and reliable—a smart investment**

EVO Telemedicine carts and wall arms are built from aircraft-grade metal and backed by a 10-year warranty.

## **Ergonomic and easy to use**

Exclusive Comfort Glide Lift Technology™ enables instant, silent, one-handed height adjustment. Research-based ergonomic features prevent unhealthy postures, awkward movements and repetitive motion injuries (RMIs).

# Popular Telemedicine Cart and Wall Mount Configurations

All EVO Telemedicine workstations are fully customizable to meet the requirements of any telehealth workflow



## EVO LCD Telemedicine Cart

- EVO LCD cart with 500 Wh LiFePO<sub>4</sub> power system
- VESA camera mount
- Speaker/sound bar mount
- 5" diameter wipes holder



## EVO Lite Tablet Telemedicine Cart

- EVO Lite Tablet cart (non-powered)
- Speaker/sound bar mount
- Adjustable tablet mount (locking or non-locking)
- 3" EVO ClearView drawer cabinet
- 5" diameter wipes holder



## EVO Telemedicine Cart with ClearView™ Storage

- EVO LCD cart with 500 Wh LiFePO<sub>4</sub> power system
- VESA camera mount
- Speaker/sound bar mount
- 3"/3"/3" EVO ClearView drawer cabinet (eLock)
- 1 large bin and 3 mini bins
- 5" diameter wipes holder



## EVO Telemedicine Cart with Bin Storage

- EVO LCD cart 500 Wh LiFePO<sub>4</sub> power system
- VESA camera mount
- Speaker/sound bar mount
- 3"/3"/3" med bin storage cabinet (e-Lock)
- Standard storage bin kits (2)
- Wide storage bin kit (1)
- 5" diameter wipes holder



### **EVO Telemedicine Wall Mount**

- EVO WA-JT wall arm workstation
- VESA camera mount
- Speaker/sound bar mount
- Digital signature pad holder
- Barcode scanner mount



### **EVO Telemedicine Wall Cabinet**

- EVO WS-27 wall cabinet
- VESA camera mount
- Speaker/sound bar mount




# Customize to create exactly the solution you need

Add EVO accessories to support a full range of telehealth workflows

## ClearView Storage Options



**Single Drawer**  
14.5"w x 9.25"d x 4.5"h



**Double Drawer**  
14.5"w x 9.25"d x 8.25"h



**Triple Drawer**  
14.5"w x 9.25"d x 11.9"h



**Single Drawer**  
12.5"w x 7.5"d x 3.5"h



**1 Large Bin**  
10"w x 7"d x 3.2"h



**2 Long Bins**  
10.5"w x 3"d x 3"h



**3 Mini Bins**  
3.5"w x 7.3"d x 2.75"h



**Deep Drawer**  
12.5"w x 7.5"d x 7.1"h

## Med Bin Storage Options



**Triple-level Cabinet**  
8.14"w x 13.25"d x 12.85"h




**Standard Med Bins**  
14.5"w x 9.25"d x 8.25"h



**Wide Med Bins**  
6.75"w x 10.5"d x 3"h

\* All drawer and bin dimensions are interior dimensions.



Camera mounts



Speaker/soundbar mounts



Label/wristband printer mounts



Barcode scanner mounts




Digital signature pad holder



Glove box holders



Baskets and wipes holders



Worksurface utility cup

## Add videoconferencing equipment from leading manufacturers (partial list):

Bose | Cyber Acoustics | Dell | Jabra | JBL | Lenovo | Logitech | Lumens Digital Optics | Panasonic | Poly (a division of HP) | PTZ Optics | Sony | Vaddio

## Collaborate with us to design custom mounts for:

Medical and dental cameras | Speakerphones | Stethoscopes | Vital sign monitors | Spirometers | ECG monitors | Ultrasound probes | Oximeters | Scopes

## Create custom telehealth workstations for:

Remote Patient Screening | Quarantine Care Engagement | Interpretive Services | Telestroke Treatment | Telepsychiatry | Virtual Rounding  
Virtual Care | Telepresence | Specialty Consultations



NOW PART OF  
**GCTX**  
Healthcare Solutions

140 Constitution Blvd. | Franklin, MA 02038  
800.649.2278 | sales@gcx.com